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KT24 Sp. z o.o.Sp.k
Sikorskiego 86C
63-100 Śrem
NIP: 7851807321

INTRODUCTION

Read this manual before assembling, installing, and using the product. Follow it to avoid situations that may lead to damage to the device or injury or even death of the user and people in the vicinity. The manufacturer is not responsible for damage caused by improper use of the equipment.

Keep this manual for future reference.

SAFETY RULES

- Use this device only as indicated in these operating instructions. Any deviation from the standard is dangerous. The manufacturer is not responsible for damages that may result from improper use.
- Carefully unpack the device. Keep all materials for possible future reference. Do not allow children to play with the packaging materials. Check the product for damage. Do not attempt to use or repair it if any part appears damaged. Refer servicing to a qualified person.
- Before connecting, check the nameplate to see if the power is the same as the household power. Do not connect if the technical characteristics are different. The manufacturer cannot be held responsible for damage caused by improper use and failure to observe this rule.
- Before each use, check that the power cord is not damaged. During use, check that the cable does not come into contact with hot parts of the equipment. Do not use the equipment if the power cord is damaged.
- If the power cord is damaged, you must have it repaired by an authorized service center.
- If you smell a strange odor while using the equipment, turn it off immediately and unplug it to avoid damage.
- This appliance must not be used by children under 8 years of age and by persons with limited sensory abilities or lack of experience. Such persons may use the equipment only after you have instructed them on how to use it safely and they understand the hazards.

- Children must not play with the equipment.
- Cleaning and maintenance must not be carried out by children without supervision.
- Do not use the equipment with wet hands or bare feet.
- Never leave the equipment outside due to weather conditions (rain, sun, cold, etc.)
Store it in a dry place only.
- This equipment is designed for home use only and must not be used for commercial purposes.
- Always unplug the equipment and wait until it has cooled down before cleaning.
- Never use the device if it is damaged or in case of a short circuit. Do not attempt to repair the device under any circumstances. Have it repaired by an authorized service center. Use only original parts.
- Do not drop or hit the equipment.
- Do not use the equipment if it has been dropped.
- Protect the equipment from liquid. Never immerse in water or other liquids. If the equipment falls into water, unplug it immediately and contact authorized service personnel.
- Make sure that the plug is securely connected to the power outlet to avoid overheating.
- When you pull the plug, never pull on the cord.
- Maximum storage temperature for the equipment is 40 degrees Celsius.
- The unit can be operated outdoors during dry weather. It must not get wet.

RESIDUAL RISK

Even if the appliance is used for its intended purpose and all guidelines in these instructions are followed, it is not possible to eliminate residual risks.
Use common sense when operating the device. Use caution.

LIST OF COMPONENTS AND CONSTRUCTION

1. Massaging tip
2. Massage head connector
3. LCD display
4. Battery Box
5. Handle

BEFORE FIRST USE

1. Remove all packaging material.
2. Remove all remnants of packaging and shipping protection.
3. Make sure the package contents are complete.
4. Check the unit, the power cord, the power plug, and all additional parts for possible damage during shipping.
5. If possible, store the packaging components.
6. Dispose of in accordance with an accepted waste disposal system.

DEVICE OPERATION

1. Fully charge the battery up to six hours before first use.
2. To charge the battery, plug the DC end of the supplied 24V adapter into the battery charging port and plug the AC side cord into a wall outlet.
3. The LED bar on the battery flashes to show the battery charge level and to indicate active charging.
4. With the battery connected, the LEDs on the unit flash to show the battery charge level and indicate active charging.
5. A full charge is indicated when all LEDs stop flashing and remain fully illuminated.
6. The battery can be charged at any time and at any charge level.
7. It is not recommended to completely exhaust the battery to the level of the red LED.

8. Average run time is 2+ hours depending on speed level and pressure level during use.
9. Align and insert the battery into the base of the unit - make sure it is securely locked in place.
10. Turn on the device by sliding the power switch on the bottom of the battery to the ON position, press the power button once to turn on the device to level three. Pressing the power button again will turn off the device.
11. With the power switch in the ON position, press the power button once to turn the unit to level one, twice to turn the unit to level two, and three times to turn the unit to level three. Pressing the power button again will turn the unit off.
12. Blue LEDs indicate the appropriate power level.
13. When the device is off, insert the applicator completely by pushing it in. Remove it by pulling it out.
14. Massage the desired body part at the appropriate speed, applying pressure as needed.
15. To remove the battery, lightly press the battery release button with your thumb. Use your other hand to firmly pull down on the battery holder.
16. With your other hand, pull down on the battery holder.

SPECIFICATION

Voltage and frequency	110-240v ~ 50/60 Hz
Battery capacity	25,2V – 2500 mAh
Speed regulation	1200-3500rpm
Noise level	55dB
Work time	4,5-5h
Dimensions	238mmx258mm
Weight	800g

CLEANING AND MAINTENANCE

1. Use a slightly damp towel to clean the unit, then wipe it dry with a soft cloth.
2. For storage or travel, remove the battery and store it in the carton provided.
3. Do not immerse in water. Keep away from liquids and heat sources. Keep ventilation ports free of dust and debris.
4. Do not remove screws or attempt to disassemble the unit.
5. Do not operate continuously for more than one hour. Let the unit stand for 30 minutes before using it again.
6. Disconnect the charger from the device battery after charging or before use.
7. For adult use only. Before using this product, consult your physician.

UTILIZATION

The equipment should not be disposed of in the same way as other waste. To prevent the harmful effects of waste on the environment and human health, it is prohibited to place used equipment together with other waste.

Each user is obliged to return such equipment to the collection point.

GUARANTEE

The terms of the guarantee applicable in each country have been published by our respective national distributors. We will repair possible faults on your unit free of charge within the warranty period, insofar as the faults are caused by material defects or defective workmanship. Guarantee claims should be addressed to your dealer or the nearest authorized customer service centre, and supported by documentary evidence of purchase.

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